

METS Coned Hemi-Pelvis

- 
-
- Evolved from the experience of more than 250 custom-made hemi-pelvic replacements
 - HA coated stem for improved bone fixation
 - Large internal diameter shell for cup sizes up to 50mm outer diameter
 - Selection of 4 stem diameters with cutting flutes for rotational stability
 - Axisymmetric hole array for optimal screw alignment and cementing
-




Indications

Major acetabular bone loss
Revision hip procedures
Bone tumours

Stanmore Implants

210 Centennial Avenue
Centennial Park
Elstree
WD6 3SJ
United Kingdom

T +44 (0) 20 8238 6500
F +44 (0) 20 8953 0617
www.stanmoreimplants.com

Coned hemi-pelvis		Stem				
	Ø9>7mm	Length=93mm	●		mchp/Ø9	
	Ø10>8mm	Length=93mm	●		mchp/Ø10	
	Ø12>9mm	Length=93mm	●		mchp/Ø12	
	Ø16>10mm	Length=93mm	●		mchp/Ø16	
Steinmann						
bone pins (x3)	Ø5mm	Length=300mm	●		msp/300Ø5	
Bone screws (x3)						
	Ø4.5mm	Length=50mm	●		mbscrwx/50Ø4.5	
	Ø4.5mm	Length=80mm	●		mbscrwx/80Ø4.5	