

# METS Modular Total Femur

- Evolved from the experience of more than 3000 femoral replacements
- SMILES femoral component specifically designed for limb salvage application
- 3 tibial hinge configurations
- Range of trochanteric components
- Precision-made tapered junctions ensuring integrity and durability

## Indications



































































Primary bone tumours  
Failed joint replacements  
Failed massive allografts  
Failed massive replacements

## Stanmore Implants

210 Centennial Avenue  
Centennial Park  
Elstree  
WD6 3SJ  
United Kingdom

T +44 (0) 20 8238 6500  
F +44 (0) 20 8953 0617  
[www.stanmoreimplants.com](http://www.stanmoreimplants.com)



<b>Trochanters</b>							
<b>Smooth</b>	Small	left		mstrc/LSmlU			
<b>Uncoated</b>	Small	right		mstrc/RSmIU			
	Standard	left		mstrc/LStdU			
	Standard	right		mstrc/RStdU			
<b>Spiked</b>	Small	left		mstrc/LSmlC			
<b>HA coated</b>	Small	right		mstrc/RSmIC			
	Standard	left		mstrc/LStdC			
	Standard	right		mstrc/RStdC			
<b>Trochanter plate</b>				msfpte			
<b>Trochanter screws</b>	Short			msfscw/Short			
	Medium			msfscw/Medium			
	Long			msfscw/Long			
<b>Trochanter reattachment wire</b>	0.6m			mstrw			
<b>Total femoral link components</b>	165mm			mstfl/165			
	225mm			mstfl/225			
<b>Principal shafts</b>	45mm			msfshft/45			
	60mm			msfshft/60			
	75mm			msfshft/75			
	90mm			msfshft/90			
	105mm			msfshft/105			
	120mm			msfshft/120			
	135mm			msfshft/135			
	150mm			msfshft/150			
<b>Femoral knees</b>	Small	Left		mkfe/LSm			
	Small	Right		mkfe/RSm			
	Standard	Left		mkfe/LStd			
	Standard	Right		mkfe/RStd			
<b>Tibial: rotating hinges</b>	Small			mkrhp/Sm			
<b>Polyethylene</b>	Standard			mkrhp/Std			
<b>Tibial: rotating hinges Metal casing</b>	Small	Short Stem		mkrhm/SmSt			
	Standard	Short Stem		mkrhm/StdSt			
	Small	Long Stem		mkrhm/SmLg			
	Standard	Long Stem		mkrhm/StdLg			
<b>Tibial: fixed hinges</b>	Small	Short Stem		mkfh/SmSt			
	Standard	Short Stem		mkfh/StdSt			
	Small	Long Stem		mkfh/SmLg			
	Standard	Long Stem		mkfh/StdLg			
<b>Tibial: plateau plates</b>							
<b>Small</b>	Sml 5mm			mktp/Sm5			
	Sml 10mm			mktp/Sm10			
	Sml 15mm			mktp/Sm15			
	Sml 20mm			mktp/Sm20			
<b>Standard</b>	Std 5mm			mktp/Std5			
	Std 10mm			mktp/Std10			
	Std 15mm			mktp/Std15			
	Std 20mm			mktp/Std20			
<b>CoCr femoral heads</b>	Ø28mm	-3.5mm		msfmh/cc28-3.5			
	Ø28mm	0mm		msfmh/cc28-0			
	Ø28mm	+3.5mm		msfmh/cc28+3.5			
	Ø28mm	+7mm		msfmh/cc28+7			
	Ø28mm	+10.5mm		msfmh/cc28+10.5			
	Ø32mm	-4mm		msfmh/cc32-4			
	Ø32mm	0mm		msfmh/cc32-0			
	Ø32mm	+4mm		msfmh/cc32+4			
	Ø32mm	+8mm		msfmh/cc32+8			
<b>Ceramic femoral heads</b>							
<b>Large femoral heads Ø34-56mm in 2mm increments</b>					<b>Available on special request</b>		